

Test Menu

Category	Abbr.	Assays
Diabetes	HbA1c	Glycosylated Hemoglobin
Thyroid	TSH	Thyroid Stimulating Hormone
	TT3	Total Triiodothyronine
	TT4	Total Thyroxine
	FT4	Free Thyroxine
	☆ FT3	Free Triiodothyronine
Fertility	FSH	Follicle-Stimulating Hormone
	PROG	Progesterone
	LH	Luteinizing Hormone
	AMH	Anti-Mullerian Hormone
	β-HCG	β-Human Chorionic Gonadotropin
	☆ E2	Estradiol
	☆ PRL	Prolactin
	☆ T	Testosterone
	IL-6	Interleukin-6
Infection	CRP	C- Reactive Protein
	HBP	Heparin Binding Protein
	PCT	Procalcitonin
	S100-β	Central Nervous System Specific Protein
Craniocerebral	Hs-cTnI	High Sensitivity Cardiac Troponin I
	D-Dimer	D-Dimer
	ST2	Soluble Growth Stimulating Expression Gene 2
Cardiovascular	25-OH-VD	25-hydroxy-Vitamin D
	VB12	Vitamin B12
Microelement	Folic Acid	Folic Acid
	☆ Fer	Ferritin

☆ Coming soon

R&D Center: 3938 Trust Way, Hayward, CA 94545, USA

Production Base: No.3 Building, Fangyi Science Park,
No.365 Huai An Road, Shijiazhuang, Hebei, China.

Website: www.hipro.us

Email: hipro@hipro.us



PalmF

Fluorescence Immunoassay Analyzer

Portable but Powerful
Multi-scenario application

Rapid Detection
Test in 155 seconds

Long-lasting Endurance
Up to 584 consecutive tests

Hipro

Hipro Biotechnology Co., Ltd.



High Precision

30 times
of sensitivity



Wide Linearity

3 times of industry average

High Sensitivity

Unique optical path design

Dual Calibration


Optical path & curve calibration

Long-lasting Power

584 consecutive tests for
AAA batteries



Multi-scenario Application



Your powerful laboratory in hand

Smart System

Intelligent interpretation

Avoid manual interpretation errors



Acousto-optic system

Avoid operation errors
Two-color indicator light
Smart power management



One-step Operation

Add Sample



Reaction



Read the result



Portability and Flexibility

Specimen types include serum, plasma and whole blood

Exquisite design, suitable for one-handed operation

Parameters

Name	PalmF Fluorescence Immunoassay Analyzer
Size	172 × 80 × 56mm
Weight	≤ 220g (exclude battery)
Power	≤1.5VA